


STANDARD OPERATING PROCEDURE OF DC POWER SUPPLY

 <p>BLDEA'S Shri Sanganabasava Mahaswamiji College of Pharmacy & Research Centre Vijayapur 586103</p>	DEPARTMENT : PHARMACEUTICS	SOP NO: Company Name: servo controlled voltage stabilizer Model no: ARPS 302D
	AREA: PHARMACEUTICS LAB	PAGE NO:1-2
	SUBJECT: SOP FOR DC POWER SUPPLY	EFFECTIVE DATE: 25/2/2022
		REVIEW PERIOD: 25/2/2023

Purpose:

To provide a procedure for the operating of DC power Supply apparatus

Scope:

Applicable to operation of DC Power Supply apparatus

References:

Instrument Manual

Safety Issues and Precautions:

- Were Head cap, Mask, Hand gloves and other safety requirements during the performance of stated activity
- Please make the entry of Usage with require details in the Instrument log book
- Handle the Instrument properly with care
- Always connect a ground wire.
- Not touch high voltage areas.
- Cover high voltage areas.
- Share an awareness of danger.
- Perform operations with your right hand.
- Turn off the power before touching equipment.
- Pay attention to electric charge in cables.

Procedure:

- 1) Put on all safety equipment and tie back loose hair or clothing before operating.
- 2) Check the make sure all parts of machine are working and functioning properly.
- 3) Make sure power supply is off and set up wires and complete circuit as desired.
- 4) When ready, turn ON power supply.
- 5) When finished, turn OFF power supply before disconnecting any wires.
- 6) Disconnect circuits and place all wires and other equipment back in the appropriate place.
- 7) Clean workspace.

	PREPARED BY	CHECKED BY	APPROVED BY
NAME	Mrs. Ashwini S. Gudigenavar	Mr. A B. Walikar	Dr. S C. Marapur
DESIGNATION	Asst. professor	Asst. Professor	Professor & Head
DEPARTMENT	Pharmaceutics	Pharmaceutics	Pharmaceutics
SIGNATURE & DATE			