

 <p>BLDEA'S Shri Sanganabasava Mahaswamiji College of Pharmacy & Research Centre Vijayapur 586103</p>	DEPT: PHARMACEUTICAL CHEMISTRY	SOP NO.: BCP/PC/SOP/017
	INSTRUMENT: MELTING POINT APPARATUS	PAGE NO.: 01-02
	MAKE: ConTech MODEL: CDMP-300-01 PROCURED ON: 10-08-2020	EFFECTIVE DATE: 01/01/2022
	SUBJECT: SOP FOR MELTING POINT APPARATUS	REVIEW PERIOD: 31/12/2022

Objective:

The following document describes the standard operating procedure for Melting point apparatus.

Scope:

Melting Point Apparatus measures the melting point of a solid substance i.e. the temperature at which the state of the substance changes from solid to liquid. This is an indication of the purity of a substance

Procedure:

1. Plug in to ensure the power supply.
2. Switch **ON** the main power supply and instrument mains.
3. The screen displays the starting temperature and rate of heating as per the previous setting (setting before the instrument was switched off last time).
4. Initial temperature can be changed by entering required temperature through the keyboard and the press **ENTER** key (↵).
5. Cursor moves to “**RATE**” (heating rate). Input the rate of heating in °C per minute and the press **ENTER KEY** (↵).
6. The temperature starts rising at a higher heating rate till the set value of starting temperature & settles at the initial temperature value.
7. At this stage insert the capillary filled with sample upto 3mm height (for consistency of melting point results loading height should be consistently maintained) & press **START** key.
8. Now the temperature starts rising slowly as per the set value of rate of heating.

9. The screen displays the graph of light transmittance(y axis) v/s melting point temperature (x axis) with initial melting point (S Melt) & Final melting point (F Melt) in °C. The current value of temperature at capillary is also displayed on the top of the screen with ↑indicating increase in the temperature. The screen records the initial and final melting values along with the graph. After the detection of final melting value the heating stops & the temperature settles at the set initial temperature.
10. In case the capillary is not inserted while starting the test (pressing **START** key) the screen displays an error.
11. This instrument allows maximum 3 tests to be carried for finding average value of melting point results.

Precautions:

- Ensure that the capillary tube is **SEALED** on one end and **OPEN** on the other end.
- **DO NOT** fill the capillary tube more than 25% of its capacity.

	PREPARED BY	CHECKED BY	APPROVED BY
NAME and DESIGNATION	Ms.Hasti.Kenia Assistant professor	Dr. S.M.Metri Associate Professor	Dr.B.Shivakumar Professor & Head
SIGNATURE & DATE			