



BLDEA's SSM College of Pharmacy and Research Centre
Vijayapur - 586103
Dept. of Pharmacology

STANDARD OPERATING PROCEDURE
Updated: 01-02-2022 Valid upto 01-08-2022

TISSUEHOMOGENIZER

FUNCTION: To homogenate the tissue preparation.

Standard Operating Procedure:

- Ensure the homogenizer probe is clean before and after use.
- Add frozen tissue to homogenizer tube prior to homogenizing.
- Homogenize tissue.
- Upon completion, rinse probe with clean water.
- **Prepared by: Shri Shridharkumar.S. Biradar**
- **Checked by: Dr. V.P.Patil**
- **Approved by: Dr. Shivakumar Hugar**